



## Replacement Lamps For Dental Equipment

Consumer Price Sheet- January 1, 2025 (Good through June 30, 2025)

All Prices Are Subject to Change Without Notice

**Operatory, Curing, Fiber Optic, Handpiece, Video**

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price	Fig.
AL1824	24W Solarc Lamp	BW.WA09500 Welch Allyn Brand	---	\$541.00**	8
AL1824	24W Solarc Lamp	BW.09500-USH Ushio Brand	---	415.05**	8
AL1824	24W Solarc Lamp	BW.09500 Import Brand	---	269.00**	8
BHB	120V,250W	BW.BHB	---	44.50**	7
BU7012-GEN1	MK DENT QC5016 Coupler	BW.BU7012-GEN1 1 Pack Generic	---	47.60	42
BU7012-GEN3	MK DENT QC5016 Coupler	BW.BU7012-GEN3 3 Pack Generic	---	128.95	42
DNE	120V,150W	BW.DNE	3019	41.10	6
DNF	21V,150W	BW.DNF	3023	41.10	6
DZE	24V,150W	BW.DZE Osram Brand	---	30.35**	17
EFM	8V,50W	BW.EFM	3040	31.45**	2
EFN	12V,75W	BW.EFN	3035	28.99	1
EFN OSRAM	12V,75W	BW.EFN-OSR Osram Brand	3035	31.98	1
EFP	12V,100W	BW.EFP	3020	28.70	1
EFP OSRAM	12V,100W	BW.EFP-OSR Osram Brand	3020	31.45	2
EFR	15V,150W	BW.EFR	3014	28.70	1
EFR OSRAM	15V,150W	BW.EFR-OSR Osram Brand	3014	31.50	1
EHJ	24V,250W	BW.EHJ	3038	15.07	14
EHT	120V,250W	BW.EHT	4073	23.25**	33
EJA	21V,150W	BW.EJA	3018	22.76	1
EJA OSRAM	21V,150W	BW.EJA-OSR Osram Brand	3018	30.28	1
EJA PHILIPS	21V,150W	BW.EJA-PHIL Philips Brand	3018	30.84	1
EJL	24V,200W	BW.EJL	3017	27.85	1
EJL OSRAM	24V,200W	BW.EJL-OSR Osram Brand	3017	29.25	1
EJL USHIO	24V,200W	BW.EJL-USH Ushio Brand	3017	29.92	1
EJM	21V,150W	BW.EJM	3011	28.54	1
EJM OSRAM	21V,150W	BW.EJM-OSR Osram Brand	3011	31.11**	1
EJV	21V,150W	BW.EJV	3022	28.01	1
EJV EIKO	21V,150W	BW.EJV-EIKO Eiko Brand	3022	9.00**	1
EJV OSRAM	21V,150W	BW.EJV-OSR Osram Brand	3022	30.17	1
EJV USHIO	21V,150W	BW.EJV-USH Ushio Brand	3022	33.37	1
EKE	21V,150W	BW.EKE	3027	19.15	1
EKE GE	21V,150W	BW.EKE-GE GE Brand – No Longer Made By GE – See Other Brands			
EKE USHIO	21V,150W	BW.EKE-USH Ushio Brand	3027	23.84	1
ELC	24V,250W	BW.ELC	---	23.27	2
ELZ	21V,150W	BW.ELZ	3024	40.34**	6
ENL	12V,50W	BW.ENL	3036	21.37**	1
ENZ	30V,50W	BW.ENZ	4084	35.60**	1
EPN	12V,35W	BW.EPN	---	44.34	1
EPV	14.5V,90W	BW.EPV	4080	36.20**	1
EPX	14.5V,90W	BW.EPX	---	34.50**	1
ESD	120V,150W	BW.ESD	3028	38.52	1
ESD USHIO	120V,150W	BW.ESD-USH Ushio Brand	3028	49.95	1
ESJ	82V,85W	BW.ESJ	---	39.25**	1
ESP	120V,150W	BW.ESP	---	27.20**	19
ESX	12V,20W	BW.ESX	4085	13.64	1
ESY	100V,150W	BW.ESY	3034	30.15**	13
ETJ	120V,250W	BW.ETJ	3015	31.77	1
ETJ OSRAM	120V,250W	BW.ETJ-OSR Osram Brand	3015	34.45	1
EVA	12V,100W	BW.EVA	3059	18.64	12
EVA OSRAM	12V,100W	BW.EVA-OSR Osram Brand	3059	21.73	12
EZK	120V,150W	BW.JCR120V150WB	3008	40.53	1

\*\*Price Is For Stock On Hand Only Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!  
E-mail: [bulbwork@bulbworks.com](mailto:bulbwork@bulbworks.com) • Order Online: [www.bulbworks.com](http://www.bulbworks.com)

<b>FCS</b>	24V,150W	<b>BW.FCS</b>	4070	\$10.77	14
<b>FCS OSRAM</b>	24V,150W	<b>BW.FCS-OSR</b> Osram Brand	4070	11.99	14
<b>FDS/DZE</b>	24V,150W	<b>BW.FDS.DZE</b>	4074	22.60	17
<b>FDS/DZE</b>	24V,150W	<b>BW.FDS.DZE-OSR</b>	4074	24.70	17
<b>FDS/DZE</b>	24V,150W	<b>BW.FDS.DZE-USH</b>	4074	25.24	17
<b>FDT</b>	12V,100W	<b>BW.FDT</b>	4081	18.77	18
<b>FTD</b>	12V,20W	<b>BW.FTD</b>	---	12.04	5
<b>JA12V-55WD H3</b>	12V,55W Flat Connector	<b>BW.JA12V-55WD Flat Connector</b>	4071	13.91	28
<b>JAV12V-55W</b>	12V,55W Pin Connector	<b>BW.JA12V-55W Pin Connector</b>	-	15.75	-
<b>JA12V-100WD</b>	12V,100W	<b>BW.JA12V-100WD</b>	---	13.74**	28
<b>JA24V-60WDBT</b>	24V,60W	<b>BW.JA24V60WDBT</b>	---	35.33	29
<b>JA24V-70WD</b>	24V,70W	<b>BW.JA24V-70WD</b>	---	18.44	28
<b>JC6V-20W</b>	6V,20W	<b>BW.JC6V-20W</b>	---	17.87	11
<b>JC12V-50W</b>	12V,50W	<b>BW.JC12V-50W</b>	3032	14.11	12
<b>JC12V-50WG1</b>	12V,50W	<b>BW.JC12V-50WG1</b>	---	Call	12
<b>JCD24V-60WBA</b>	24V,60W	<b>BW.JCD24V-60WBA</b>	---	74.37	19
<b>JCD24V-75W</b>	24V,75W	<b>BW.JCD24V-75W</b>	---	35.87	19
<b>JCD120V-150WB</b>	120V,150W	<b>BW.JCD120V150WB</b>	---	36.20**	12
<b>JCR12V75WMR8</b>	12V,75W	<b>BW.JCR12V75WMR8</b>	---	66.90**	-
<b>JCR120V-150WB</b>	120V,150W	<b>BW.JCR120V150WB</b>	3008	40.53	1
<b>JCRM12V-50W</b>	12V,50W	<b>BW.JCRM12V-50W</b>	---	35.90**	3
<b>JCRM12V-75W</b>	12V,75W	<b>BW.JCRM12V-75W</b>	3016	37.59	4
<b>JCRM12V-100W</b>	12V,100W	<b>BW.JCRM12V-100W</b>	3021	41.37	4
<b>M10P003</b>	14W/9.5W MR11 Solarc Lamp	<b>BW.M10P003</b>	---	155.94	7
<b>M21E001</b>	60V,21W MR11 Welch Allyn Brand	<b>BW.M21E001</b>	---	Call	9
<b>M21E001-IMP</b>	60V,21W MR11 Import Brand	<b>BW.M21E001-IMP</b>	---	235.29**	9
<b>M34 (Delisted)</b>	6V,20W	<b>BW.JC6V-20W</b>	---	17.87	11
<b>MHR150N</b>	95V,150W	<b>BW.MHR150N--Discontinued</b>	---	Call	10
<b>PELTON &amp; CRANE LAMP TOOL</b>	To Install Lamp	<b>BW.PC TOOL</b>	4078	12.90	31
<b>Q75CL/DC-28V</b>	28V,75W	<b>BW.Q75CL.DC-28V</b>	4079	26.60	19
<b>Q150T4/CL25V</b>	25V,150W	<b>BW.Q150T4.CL25V</b> Ushio Brand	4076	24.66	30
<b>Q150T4/CL25V IMP</b>	25V,150W	<b>BW.Q150T4CL25IMP</b> Import Brand	4076	20.10	30
<b>SABLE 1</b>	Access L & Rotomax L Handpiece	<b>BW.SABLE-1</b>	---	Call	36
<b>SABLE 5</b>	Access L & Rotomax L Handpiece	<b>BW.SABLE-5</b>	---	Call	36
<b>THE71</b>	15V,150W	<b>BW.EFR</b>	---	28.70	2
<b>VRB-GEN1</b>	Vector VKS6 Coupler Bulb	<b>BW.VRB-GEN1</b> 1 Pack Generic	---	Call	41
<b>VRB-GEN3</b>	Vector VKS6 Coupler Bulb	<b>BW.VRB-GEN3</b> 3 Pack Generic	---	Call	41
<b>VSBL</b>	LED For Sirona 6 Pin Coupler	<b>BW.VSBL</b>	---	75.91	42A
<b>WA09500</b>	24W Solarc Lamp	<b>BW.WA09500</b> Welch Allyn Brand	---	541.00**	8
<b>WA09500</b>	24W Solarc Lamp	<b>BW.09500-USH</b> Ushio Brand	---	415.05**	8
<b>WA09500</b>	24W Solarc Lamp	<b>BW.09500</b> Import Brand	---	269.00**	8
<b>WA09800</b>	21W Metal Halide	<b>BW.WA09800</b>	---	Call	8
<b>Y900529</b>	NSK Handpiece Bulb 1 Pack	<b>BW.Y900529-GEN1</b> Generic Brand 1 Pk	---	54.54**	41
<b>Y900529</b>	NSK Handpiece Bulb 3 Pack	<b>BW.Y900529-GEN3</b> Generic 3 Pk	---	126.27**	41
<b>3S6-120V</b>	120V,3W	<b>BW.3S6-120V</b>	---	4.50	21
<b>25A-120V</b>	120V,25W	<b>BW.25A-120V</b>	---	6.85	22A
<b>100A21SP</b>	120V,100W	<b>BW.100A21SP</b>	---	Call	22B
<b>125T10/P</b>	120V,125W	<b>BW.125T10.P</b>	4072	71.38	20
<b>300SB/IF</b>	120V,300W	<b>BW.300SB.IF</b>	---	30.72**	26
<b>1222 LED</b>	Kinetic Handpiece - LED Version - Rated 3.5V	<b>BW.1222-LED</b> Can Be Used To Replace Standard 1222 Bulb	---	78.50	44
<b>13165</b>	14V,35W	<b>BW.13165</b> Philips Brand - No Longer Made By Philips - See BW.13165-IMP	---	---	---
<b>13165 IMPORT</b>	14V,35W	<b>BW.13165-IMP</b> Import Brand	3026	33.10	3
<b>13298</b>	10V,52W	<b>BW.3367-5267**</b>	3009	24.08	3
<b>13865</b>	12V,75W	<b>BW.13865</b>	---	47.89**	4
<b>14552</b>	12V,75W	<b>BW.14552</b>	---	6.00**	3
<b>21237</b>	12V,80W	<b>BW.21237</b>	---	45.63	4
<b>22076 LED</b>	W & H Handpiece - LED Version - Rated 3.5V	<b>BW.22076-LED</b> Can Be Used To Replace Standard 22076 Bulb	---	77.87**	45

60986	2 Pack Star Quick Connect Swivel	BW.60986	---	\$68.50**	32
64115	12V,10W	BW.64115	---	16.70**	55
64156	24V,70W	BW.64156	---	19.67	28
64440	12V,50W	BW.64440	---	10.76	12
64450AX	12V,75W	BW.64450AX	---	13.59	15
64451	24V,60W No Black Tip	BW.64451	---	37.70	28
64602	12V,50W	BW.64602	---	32.75	12
64613	12V,75W	BW.64613-IMP	3010	33.01	4
64617	12V,75W	BW.64617	3016	38.35	4
64617S	12V,75W	BW.64617S	---	45.60	4
64624	12V,100W	BW.64624	3021	41.95	4
64638	24V,100W	BW.64638	---	29.71	14
64638	24V,100W	BW.64638-IMP	---	21.52	14
95775	41V,200W	BW.95775	---	69.42	1
262522	3.3V Halogen Star 6 Pin Swivel Coupler	BW.262522	---	41.10	44A
265385	3.3V LED Star 6 Pin Swivel Coupler	BW.265385	---	78.25	48A
561398	13V,100W	BW.561398	---	72.05**	3
650064	2 Pk Midwest Power Optic For In-Hose Swivel	BW.650064-GEN2 Generic Brand 2 Pk	---	79.10**	38
650181	DCI,Dentsply,Kinetic,Midwest,Schein Light-Pak	BW.650181	---	122.00**	43A
650197	3 Pk Midwest Power Optic Handpiece Bulb	BW.650197 Midwest Brand	3025	Call	35
790257	Midwest Stylus Handpiece	BW.790257GEN1	---	27.95**	36
790257 LED	Midwest Stylus Handpiece – LED Version	BW.790257-LED Can Be Used To Replace Standard 790257 Bulb	---	97.34	46
553-3881	2 Pk Kavo Lux Coupler Bulb – 3.3V,4W	BW.553-3881 Kavo Brand – 2 Pack	---	171.75**	42
553-3881	1 Pack Kavo Lux Coupler Bulb – 3.3V,4W	BW.5533881-GEN1 Generic Brand – 1 Pack	---	47.60	42
553-3881	3 Pack Kavo Lux Coupler Bulb – 3.3V,4W	BW.5533881-GEN3 Generic Brand – 3 Pack	---	128.95	42
553-3881 LED	Kavo Lux Coupler Bulb – LED Version - Rated 3.3V, 4W	BW.553-3881-LED Can Be Used To Replace Standard 553-3881 Bulb	---	81.90	47
0.01905510 LED	ADEC LT25 Low Speed Motor	BW.0.01905510LE	---	92.95	48
0.02207610 LED	ADEC High Speed Handpiece	BW.0.02207610LE	---	Call	45A
0.0352301	ADEC RA24 Roto-Quick Coupler	BW.0.03523010	---	92.95	43
0.04300400	ADEC Handpiece	BW.0.04300400	---	99.80	47A
041-117-00	ADEC Excellence (Older Style)	BW.041-117-00	---	Call	39
041-317-00	ADEC "Newer" Style 3.5V	BW.041-317-00	3031	40.38	39
041-317-00 LED	ADEC "Newer" Style – LED Version – Rated 3.5V	BW.041-31700LED Can Be Used To Replace Standard 041-317-00 Bulb	---	97.33	48
041-709-00	17V,95W	BW.041-709-00	---	31.60	16
041-709-00	17V,95W	BW.041-709-00IM	---	23.74	16
3367-5267	10V,52W	BW.3367-5267**	3009	24.08	3
42001-02121	Dabi Atlante Generic Brand 1 Pack	BW.42001021GEN1	---	9.00**	42
42001-02121	Dabi Atlante Generic Brand 3 Pack	BW.42001021GEN3	---	25.00**	42
70201102004	13V,75W	BW.70201102004 DELISTED	---	Call	2

### Fluorescents For X-Ray Viewing; Magnifying Lights, Curing, Work Lights

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price	Fig.
CF9DS.E.841	9W Comp. Fluor. 4100K 4 Pin Base	BW.CF9DS.E.841	---	\$20.87	60
CF9DS.BLUE	9W Blue Comp. Fluor.	BW.CF9DS.BLUE	---	31.03	59
DS9W.78	9W UV Compact Fluor.	BW.DS9W.78	---	22.15	59
F6T5/CW	6W 9" Cool White 5/8" Diam.	BW.F6T5.CW	3050	9.80	61
F6T5/D	6W 9" Daylight 5/8" Diam.	BW.F6T5.D	3048	9.90	61
F8T5/CW	8W 12" Cool White 5/8" Diam.	BW.F8T5.CW	3051	10.82	61
F8T5/D	8W 12" Daylight 5/8" Diam.	BW.F8T5.D	3049	11.12	61
F10T4/CW	10W 12 1/8" Cool White	BW.F10T4.CW	---	31.83	61
F13T5/CW	13W 21" Cool White 5/8" Diam.	BW.F13T5.CW	3052	13.67	61
F13T5/D	13W 21" Daylight 5/8" Diam.	BW.F13T5.D	---	14.41	61
F14T8/CW	14W 15" CW 1" Diam.	BW.F14T8.CW	3055	12.75	62
F14T8/D	14W 15" Daylight 1" Diam.	BW.F14T8.D	3056	14.92	62
F14T12/CW	14W 15" CW 1 1/2" Diam.	BW.F14T12.CW	3053	19.64	63
F14T12/D	14W 15" Daylight 1 1/2" Diam.	BW.F14T12.D	3054	20.73	63
F15T8/CW	15W 18" Cool White 1" Diam.	BW.F15T8.CW	3057	15.52	62
F15T8/D	15W 18" Daylight 1" Diam.	BW.F15T8.D	3058	16.07	62
FC8T9/CW	22W 8" Diameter Ring	BW.FC8T9.CW	8023	20.39	58

### X-Ray Duplicator Lamp

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price	Fig.
F8T5/BLB	8W 12" Blacklight Blue 5/8" Diam.	BW.F8T5.BLB	3049B	\$22.87	64

\*\*Price Is For Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!  
Stock On Hand Only E-mail: [bulbwork@bulbworks.com](mailto:bulbwork@bulbworks.com) • Order Online: [www.bulbworks.com](http://www.bulbworks.com)

## X-Ray Darkroom Lamps; Safelights

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price	Fig.
7C7-ORANGE	120V, 7W Orange, Candelabra	BW.7C7-ORANGE	---	\$7.08	25
7C7-RED	120V, 7W Red, Candelabra	BW.7C7-RED	---	7.08	25
7.5S11/CR	7½W Red, Med. Screw	BW.7.5S11.CR	---	8.44	24
15A/R	120-130V, 15W Red Medium Screw Base	BW.15A.R	---	18.69	23
15A15-120V	120V, 15W	BW.15A15-120V	---	5.75	22C
25A/R	120V, 25W	BW.25A.R	---	10.95	23
F15T8/R	15W 18" Red Fluor. 1" Diam.	BW.F15T8.R	---	32.55	65

### Indicator, Panel, Communicator, And Miniature Lamps (Price Each - Packed 10 To A Box)

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price	Fig.
40	6.3V, .15A T3 1/4 Miniature Screw	BW.40	---	\$3.41	-
47	6.3V, .15A Miniature Bayonet	BW.47	---	3.25	54
73	14.0V, .08A T 1 3/4 Wedge	BW.73	---	3.65	40
85	28V, .04A Wedge Base	BW.85	---	4.30	40
88	6.8V, 1.91A S8 DC Bayonet	BW.88	---	4.75	57
88 GE	6.8V, 1.91A S8 DC Bayonet	BW.88-GE GE Brand	---	8.30	57
120PSB5	120V, .025A Slide Base	BW.120PSB5	---	4.07	49
336	14V, .08A Midget Groove	BW.336	---	3.25	52
376	28V, .06A T1 3/4 Midget Flange	BW.376	---	3.50	53
381	6.3V, .20A T1 3/4 Midget Flange	BW.381	---	3.86	53
387	28V, .04A Midget Flange	BW.387	---	3.60	53
387 GE	28V, .04A Midget Flange	BW.387-GE GE Brand	---	5.82**	53
787	6V, 10W G4 Base	BW.787	---	42.00**	34
789	12V, 14W T2 3/4 G4 Base	BW.789	---	42.00**	34
1073	12.8V, 1.80A SC Bayonet	BW.1073	---	3.25	56
1819	28V, .04A T3 1/4 Miniature Bay	BW.1819	---	5.25	-
7382	14V, .08A T1 3/4 Bi-Pin	BW.7382	---	4.15	37
462420	24V, 20MA Slide Lamp	BW.462420	---	4.65	50
552450	24V, 50MA Slide Lamp	BW.552450	---	5.86	51
HPR40	6V, .67A T3 SC Mini Flange	BW.HPR40	---	3.67	55
HPR53	4V, 3.4W T3 SC Mini Flange	BW.HPR53	---	6.15	55

### Full Spectrum Lighting/Color Corrected

Full Spectrum Color Corrected Fluorescent Lamps • For The General Lighting Of Dental Offices And Laboratories

Benefits Of Full Spectrum Color Corrected Lamps Include:

- Accurate Shade And Color Matching! • Diagnostic Precision! • Visual Accuracy! • Greater Efficiency Of Personnel! • Increased Eye Comfort! • Office Area Is Made More Pleasant And Invigorating! • Helps People To Work And Feel Better!

Bulb	Description	Bulbworks Part Number	Retail Price	Fig.
F15T8/DSGN50	15W 18" Color Corrected 1" Diam.	BW.F15T8.DSGN50	\$25.68**	62
F32T8/850	32W 48" Length 1" Diameter Color Corrected	BW.F32T8.850 25 Pieces Per Case. Brand Can Vary	5.95	62
FO32T8/850	32W 48" Length, 1" Diameter Color Corrected	BW.FO32.850 30 Pieces Per Case, Sylvania or SLI Brand	11.19**	62
F40T12.C50 Replaces Standard F40T12CW	48" Length, 1½" Diameter	BW.F40T12.C50 30 Pieces Per Case	15.40**	63

Call For Any Size Or Type Of Fluorescent Lamp Not Listed. Energy Saving Versions Are Also Available.

Color Corrected Fluorescents Fit Most Standard Fluorescent Light Fixtures With No Modification Necessary.